



Where to use :	With all sampling valves and W9 hose piece.	
How to use:	Manually	
Use for:	500ml Aseptic sampling bottle	
Accessories:		
Connection hose:	0.5m PTFE (included)	
Filter:	Sterile (included)	
Material:		
Steel Parts:	AISI 316L(1.4404)	
O-rings:	Silicone	
Gasket:	PTFE and Silicone	
Key measurements		
Drawing part A:	14 mm.	0.55 In.
Drawing part B:	9 mm.	0.35 In.
Bottle thread	ISO GL45	
Net weight:	0.65 kg./pc	1.43 lb
Temp. range:	1 ⁰ C - 130 ⁰ C	34°F - 266°F
Inlet pressure range:	0 – 6 bar(g)	0 – 87 psi(g)
Bottle pressure range:	0-1,2 bar(g)	0-18 psi(g)